

KESSLER LAB-PROTEOMICS PROTOCOLS

Target Discovery Institute (TDI)
NDM Research Building
Oxford University
Old Road Campus
Roosevelt Drive
Headington
OXFORD, OX3 7FZ, UK
Tel: +44 (0) 1865 612 935
Email: Benedikt.Kessler@ndm.ox.ac.uk

IN-GEL DIGESTION

Guidelines for sample preparation

(How to protect your samples from contamination with keratin)

- TRY TO AVOID ANY CONTACT OF SAMPLES AND SOLUTIONS WITH DUST, SKIN OR HAIR
- CLEAN YOUR BENCH
- WEAR GLOVES AT ALL TIMES
- ALL REAGENTS SHOULD BE PREPARED FRESH
- USE ULTRA PURE WATER FOR ALL SOLUTIONS

SAMPLE SUBMISSION

- PROVIDE SAMPLES IN 1.5 ML TUBES (NO SCREW CAPS)
- LABEL YOUR TUBE WITH YOUR NAME, SAMPLE NUMBER AND DATE
- PLEASE FILL IN ONLINE SAMPLE SUBMISSION FORM TO PROVIDE US WITH MORE INFORMATION ABOUT YOUR SAMPLE (PROTEIN ORIGIN, TAXONOMY, CONCENTRATION)

REAGENTS

(All reagents should be prepared fresh)

Wash solution:

50 % methanol

5 % acetic acid

in MilliQ-H₂O

we suggest total volumes of 50 ml

100 mM ammonium bicarbonate (79 g/mol):

0.79 g NH₄HCO₃ / 100 ml MilliQ-H₂O

50 mM ammonium bicarbonate:

dilute 100 mM stock 1:1 in MilliQ-H₂O

1 M DTT (154 g/mol):

1.54 g dithiothreitol in 10 ml MilliQ-H₂O (aliquot, keep at -20 °C)

10 mM DTT

1.0 ml 100 mM NH₄HCO₃

10 µl 1 M DTT

50 mM iodoacetamide (185 g/mol):

0.0092 g iodoacetamide in 1 ml MilliQ-H₂O

Extraction buffer I:

50 % acetonitrile

5 % formic acid

in MilliQ-H₂O

Extraction buffer II:

85 % acetonitrile
5 % formic acid
in MilliQ-H₂O

Buffer A:

98 % MilliQ-H₂O
2 % acetonitrile
0.1 % formic acid

DAY ONE

1. Cut out the protein bands with a clean scalpel, cut band into smaller pieces (1 – 2 mm³) and place the gel pieces in a 1.5 ml tube.
2. Add **200 µl** of the wash solution and rinse the gel pieces overnight at room temperature while shaking

DAY TWO

1. Remove the wash solution and add another **200 µl** wash solution. Rinse the gel pieces for 2-3 hours at RT. Centrifuge shortly and remove the wash solution.
2. Add **200 µl** acetonitrile. Dehydrate the gel pieces for 5 min at RT. Remove acetonitrile from the sample. **REPEAT THIS STEP**
3. Add **30 µl** 10 mM DTT buffer and reduce for 30 minutes at RT. Remove the solution.
4. Add **30 µl** 50 mM iodoacetamide buffer and alkylate at RT for 30 minutes. Remove the solution.
5. Add **200 µl** acetonitrile. Dehydrate the gel pieces for 5 minutes at RT. They should be whitish. Remove the acetonitrile.
6. Rehydrate the gel pieces in **200 µl** 100 mM ammonium bicarbonate for 10 minutes at RT. Remove the solution.
7. Add **200 µl** acetonitrile. Dehydrate the gel pieces for 5 minutes at RT. Remove the solution.
8. Prepare the trypsin reagent: add **1 ml** of ice-cold 50 mM ammonium bicarbonate to 20 µg of trypsin [final concentration: 20 ng/µl]. The solution (10 x stock) should be kept on ice.
9. Dilute Trypsin stock 1:10 in 50mM ammonium bicarbonate (1x stock). Add **30 µl** (1 x stock) of the trypsin solution to the sample and rehydrate on ice for 10 min (DO NOT EXCEED!) with occasional gentle mixing.

KESSLER LAB-PROTEOMICS PROTOCOLS

10. Centrifuge the gel pieces shortly and remove the excess trypsin solution.
11. Add **5 µl** of 50 mM ammonium bicarbonate to the sample. Carry out the digestion overnight at 37 °C.

DAY THREE

1. Add **50 µl** of 50 mM ammonium bicarbonate buffer to the tube containing the gel pieces. Incubate for 10 min with occasional gentle vortex mixing. Collect supernatant in a 1.5 ml tube
2. Add **50 µl** of the extraction buffer 1 to the tube containing the gel pieces. Incubate for 10 minutes with occasional gentle vortex mixing. Add supernatant to the collection tube.
3. Add another **50 µl** of the extraction buffer 2 to the tube containing the gel pieces. Incubate for 10 minutes with occasional gentle vortex mixing. Add supernatant to the collection tube.
4. Dry the sample down completely in a vacuum centrifuge. Resuspend in 20 µl buffer A

Footnote: Step 3 and 4 (reduction alkylation) from day two can be performed prior to loading your sample onto your gel. Please contact us for a protocol!